ST.330.315

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This device is intended for a specialized procedure, which should be performed by professionals qualified in Dental Implants, with knowledge of the Neodent Guided Surgery technique. For optimum results, use the product knowing the appropriate techniques. Always apply them under appropriate conditions, also in an operating room atmosphere.

DESCRIPTION

Neodent Guided Surgery Sleeves are made of titanium alloy, in accordance with ASTM standard F136, and are available in five models: regular, narrow, Facility, mount Sleeve, and start drill Sleeve. They have cylindrical geometry with a thru-hole and external retentions for adhesion of the resin of the Surgical Guide. They are available according to the table below:

Model	Color	Diameter of the Stop (in mm)	Internal diameter (in mm)	Diameter of the body (in mm)	Height
Regular	Yellow	7	5	6	4
Narrow	Purple	4.6	3.5	4	
Facility	Yellow	4	2.9	3.4	
Mount	Yellow	3.1	1.5	2.5	7.5
Start Drill Sleeve	Yellow	3.5	2	2.7	6

INDICATIONS FOR USE

The Neodent Guided Surgery technique is indicated for surgical procedures on maxilla or mandible bones with a sequence of guided instruments, allowing for flapless surgery for installation of implants in edentulous areas.

APPLICATIONS

Neodent Guided Surgery Sleeve are attached to the Surgical Guide prototype. They are compatible with use of the Neodent Guided Surgery instruments and serve as a support to guide the implant when performing the Neodent Guided Surgery technique. They must be used according to the following table:

Model	Indications for use
Regular	Installation of implants of diameter: 3.5, 3.75, 4.0, 4.3 mm.
Narrow	Installation of implants of diameter: 3.5 mm, in regions with little mesio-distal space.
Facility	Installation of Facility implants.
Mount	Insertion of the guide mount.

Model	Indications for use
Start Drill Sleeve	Drilling of the 2.0 drill, in the Start Neodent Guided Surgery technique.

Neodent Guided Surgery Sleeves are compatible with the polymer resins manufactured using the additive manufacture process (prototyping).

CONTRAINDICATIONS

There are no contraindications, provided that it is used correctly for the purposes indicated.

HANDLING

During the virtual planning, the Sleeves must be selected according to the implant planned. The regular, narrow, and Facility Sleeves must respect the distance of 9 mm from the top of the Sleeve to the platform of the implant.

After fabrication of the Surgical Guide prototype, the Sleeves must be carefully attached according to the appropriate techniques, using adhesive glue. After attaching, perform a visual inspection and make sure of the absence of glue on the top and on the inside of the Sleeve.

Use the products of the Neodent Guided Surgery line according to their own Instructions for Use.

*NOTE: it is recommended to use ethyl cyanoacrylate-based glue.

SANITIZATION

The sanitization of the Sleeve/Surgical Guide assembly must follow the instructions of the manufacturer of the prototyped Surgical Guide.

PRESENTATION AND STERILIZATION

These products are indicated for single use and supplied non-sterile.

PRECAUTIONS

- Inappropriate planning can jeopardize the performance of the implant unit leading to failures of the system, as loss or fracturing of the implant.

- Be aware in cases patients that show signs of allergy or hypersensitivity to chemical components of the material. Titanium alloy Ti6Al4V.

- This product is of single use.

- Reuse of this product may cause:
 - adverse biological effects of residual products, microorganisms and / or substances resulting from previous uses and / or reprocessing;
 - changes in physical, mechanical and chemical properties of products, macro and micro structural, that can put in risk the desired functionality. The reuse of this product does not guarantee its safety and efficacy and disclaims any warranty of products.

- Do not use the product if the packaging is

broken.

- Do not use the product with the validity expired.
- Before each procedure, make sure the pieces are properly seated.
- Before each procedure check the conditions of the Neodent surgical instruments, always respecting their useful life. Replace the instruments if there is damage, markings deleted, sharpening jeopardized, deformation and wear.
- During the surgical procedure, make sure you don't exert lateral movements with the surgical instruments, avoiding the loss of the mounting of the Sleeve.
- Ensure that the parts are not swallowed or aspirated by the patient.
- Always use the sequence of Neodent products for the Neodent Guided Surgery technique. The use of instruments or Sleeves of other manufacturers does not ensure the perfect function of the Neodent Implant System and exempts any product warranty.
- It is the professional's responsibility to use the Neodent products according to the instructions for use.

ADVERSE EFFECTS

No adverse effects will occur, provided that the indications are complied with.

POST-OPERATIVE PRECAUTIONS AND MAINTENANCE

Instruct the patient as to the need of a professional medical monitoring after the surgery and to obey the guidelines regarding the precautions, hygiene and prescription of drugs. These guidelines are the responsibility of the professional in charge.

STORAGE CONDITIONS

This product should be stored in a clean and dry location, in a maximum temperature of 40°C and protected from direct sunlight.

DISPOSAL OF MATERIAL

After use, all the products and consumables used in the surgery for the installation of dental implants may endanger the health of anyone who handles them. Before discarding them into the environment, it is recommended to take a look at the current legislation and adhere to it.

DATE OF EXPIRATION

Written on the label.

Federal (USA) law restricts this device to sale by or on the order of a licensed dentist or physician.

Not all products are available in all countries. Please, contact the authorized distributor.

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