# ST.330.316

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\* Opcional

**Titamax CM**

	Neodent Guided Surgery 2.0	Neodent Guided Surgery 2.8	Neodent Guided Surgery 3.0	Piloto Neodent Guided Surgery 3.5	Neodent Guided Surgery 3.3	Piloto Neodent Guided Surgery 3.75	Piloto Neodent Guided Surgery 4.0
3.5 mm	✓	✓		✓			
3.75 mm	✓		✓			✓	
4.0 mm	✓		✓		✓		✓

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**Drive CM**

	Broca Neodent Guided Surgery Alvim 2.0	Broca Neodent Guided Surgery Alvim 3.5	Broca Piloto Neodent Guided Surgery 3.5	Broca Neodent Guided Surgery Alvim 4.3	Broca Piloto Neodent Guided Surgery 4.3
3.5 mm	✓	✓	*		
4.3 mm	✓	✓		✓	*

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\* Opcional

**Alvim CM**

	Broca Neodent Guided Surgery Alvim 2.0	Broca Neodent Guided Surgery Alvim 3.5	Broca Piloto Neodent Guided Surgery 3.5	Broca Neodent Guided Surgery Alvim 4.3	Broca Piloto Neodent Guided Surgery 4.3	Macho de Rosca Neodent Guided Surgery 3.5	Macho de Rosca Neodent Guided Surgery 4.3
3.5 mm	✓	✓	*			**	
4.3 mm	✓	✓		✓	*		**

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**Facility**

	Broca Neodent Guided Surgery Facility 2.0	Broca Neodent Guided Surgery Facility 2.0/2.4	Broca Piloto Neodent Guided Surgery Facility	Macho de Rosca Neodent Guided Surgery Facility
10 mm	✓	✓	*	**
12 mm	✓	✓	*	**
14 mm	✓	✓	*	**

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This device is intended for a specialized procedure, which should be performed by professionals qualified in Dental Implants, with knowledge of the Neodent Guided Surgery technique. For optimum results, use the product knowing the appropriate techniques. Always apply them under appropriate conditions, also in an operating room atmosphere.

**DESCRIPTION**

Neodent Guided Surgery Drills and Neodent Guided Surgery Pilot Drills are made of surgical stainless steel. On one end, they have a coupling for a Contra-Angle handpiece according to ISO 1797-1 - *Dental rotary instruments - Shank*.

- Titamax Neodent Guided Surgery drills, on their active end, have a cylindrical profile and cutting geometry compatible with CM Titamax Implants.
- Alvim Neodent Guided Surgery drills, on their active end, have a cylindrical profile and cutting geometry compatible with CM Alvim and CM Drive Implants.
- Facility Neodent Guided Surgery drills have cutting geometry compatible with Facility implants.

Neodent Guided Surgery Drills have laser markings to determine the depth of the drilling according to the surgical planning. The laser markings already take the real length of the implants into consideration. They have an additional length of 11 mm, compensating for the thickness of the soft tissue, the height of the Sleeve, and the height of the Drill Guide.

Neodent Guided Surgery Pilot Drills have an area of reduced diameter, without cutting power (serves as a guide for drilling), a sharp edge in the transition to the larger diameter and a “stop” that determines the depth.

Neodent Guided Surgery Drills and Neodent Guided Surgery Pilot Drills receive an application of a fine film of dark-colored carbon on the cutting surface. This film has the main purposes of:

- Increasing resistance to oxidation;
- Reducing the friction between the drill and the bone;
- Increasing the resistance to wear;
- Reducing the heat generated on the bone during the osteotomy.

Neodent Guided Surgery Drill Guides have a surgical stainless steel handle and titanium alloy drill guide, in accordance with ASTM standard F136. On the shaft, there is a marking indicating the drill to be used on the Sleeve and the colors of the guides indicate the diameter of the drill to be used, according to the table below:

Color	Drill to be used on the Drill Guide
Gray	2.0
Yellow	2.4
Purple	3.5
Blue	3.75
Green	4.0
Orange	4.3

**INDICATIONS FOR USE**

The Neodent Guided Surgery technique is indicated for surgical procedures on maxilla or mandible bones with a sequence of guided instruments, allowing for flapless surgery for installation of implants in edentulous areas.

**APPLICATIONS**

The Neodent Guided Surgery Drill is used for bone drilling in the preparation of the surgical socket for the installation of CM Titamax, CM Drive, CM Alvim, and Facility implants.

The Neodent Guided Surgery Pilot Drill is used for bone drilling in the preparation of the surgical socket, as a final drill, for installation of CM Titamax, CM Drive, CM Alvim, and Facility implants, to provide the seating of the cervical portion of the implant.


The Neodent Guided Surgery Drill Guide is a surgical instrument that serves as a guide for Neodent Guided Surgery Drills on the Surgical Guide prototyped during bone drilling, in the preparation of the surgical socket.

**CONTRAINDICATIONS**

These products have no contraindications, provided that they are used correctly for the purposes indicated.

**HANDLING**

For drilling with the Neodent Guided Surgery Drill, after mounting the Surgical Guide, select the Neodent Guided Surgery Drill corresponding to the Neodent Guided Surgery Drill that will be used, and in accordance with the model of Sleeve. Each Neodent Guided Surgery Drill has a specific Neodent Guided Surgery Drill Guide. Make sure to use the set that correctly corresponds. Fit the end of the Neodent Guided Surgery Drill inside the Sleeve of the Surgical Guide, holding it firmly by the handle. Fit the Neodent Guided Surgery Pilot Drill onto the contra-angle handpiece and set the surgical motor to drilling speed, as indicated for the implant selected. Bring the drill to the hole of the Drill Guide. Activate the motor and perform the drilling with continuous motions of insertion and removal, with plenty of irrigation. This irrigation can be manual or combined with the irrigation from the motor. The drilling must respect the bone quantity and quality, the implant selected, and also the prior virtual planning. The drilling must be in accordance with the length of the implant and the laser markings on the drill. The laser markings on the drill represent the drilling depth for the length of the implant.

Neodent Guided Surgery Titamax Drill      Neodent Guided Surgery Alvim Drill      Neodent Guided Surgery Facility Drill

ATTENTION: The figure is merely illustrative. It does not represent the real dimension and characteristics of the product. Neodent Guided Surgery drills have an additional length of 11 mm.

Nominal size	Real size (in mm)
7	18
8	19
9	20
10	21
11	22
11,5	22,5
12	23
13	24
14	25
15	26
16	27
17	28

For drilling with the Neodent Guided Surgery Pilot Drill, fit the Neodent Guided Surgery Pilot Drill on the contra-angle handpiece and set the surgical motor to drilling speed, as indicated for the implant selected. Bring the drill to the Sleeve. Activate the motor and perform the drilling with continuous motions of insertion and removal, with plenty of irrigation. This irrigation can be manual or combined with the irrigation from the motor. The drilling must respect the bone quantity and quality and also the prior virtual planning. The drilling must be in accordance with the implant selected. Perform the drilling until the Stop of the Pilot Drill touches the Sleeve. NOTE: drilling with the Neodent Guided Surgery Pilot Drill must be done without the Drill Guide.

Use the sequence of drills, according to the table below:

**Titamax CM EX**

	Neodent Guided Surgery 2.0	Neodent Guided Surgery 2.8	Neodent Guided Surgery 3.0
3.5 mm	✓		
3.75 mm	✓	*	
4.0 mm	✓	*	*

✓ Recommended sequence

\* Optional

**Titamax CM**

	Neodent Guided Surgery Drill 2.0	Neodent Guided Surgery Drill 2.8	Neodent Guided Surgery Drill 3.0	Neodent Guided Surgery Pilot Drill 3.5	Neodent Guided Surgery Drill 3.3	Neodent Guided Surgery Pilot Drill 3.75	Neodent Guided Surgery Pilot Drill 4.0
3.5 mm	✓	✓		✓			
3.75 mm	✓		✓			✓	
4.0 mm	✓		✓		✓		✓

✓ Recommended sequence

**Drive CM**

	Guided Surgery Alvim Drill 2.0	Guided Surgery Alvim Drill 3.5	Guided Surgery Pilot Drill 3.5	Guided Surgery Alvim Drill 4.3	Guided Surgery Pilot Drill 4.3
3.5 mm	✓	✓	*		
4.3 mm	✓	✓		✓	*

✓ Recommended sequence

\* Optional

**Alvim CM**

	Guided Surgery Alvim Drill 2.0	Guided Surgery Alvim Drill 3.5	Guided Surgery Pilot Drill 3.5	Guided Surgery Alvim Drill 4.3	Guided Surgery Pilot Drill 4.3	Alvim Bone Tap Guided Surgery 3.5	Alvim Bone Tap Guided Surgery 4.3
3.5 mm	✓	✓	*			**	
4.3 mm	✓	✓		✓	*		**

✓ Recommended sequence

\* Recommended for bone types I and II and optional for bone types III and IV

\*\* Recommended for bone types I and II

**Facility**

	Facility Guided Surgery Drill 2.0	Facility Guided Surgery Drill 2.0/2.4	Facility Guided Surgery Pilot Drill	Facility Guided Surgery Bone Tap
10 mm	✓	✓	*	**
12 mm	✓	✓	*	**
14 mm	✓	✓	*	**

✓ Recommended sequence

\* Recommended for bone types I and II and optional for bone types III and IV

\*\* Recommended for bone types I and II

Note 1: The correct sequence of drills and the rotation speed indicated for each implant must be respected, contributing to the success of the osseointegration of the implant. Note 2: Do not interrupt the rotation of the motor with the drill inside the surgical cavity, as this may make it difficult to remove it or may cause the drill to break. Note 3: Make sure that, during drilling, the Drill Guide is securely fit to its Stop inside the Sleeve.

**SANITATION**

These products must be correctly sanitized before each use.

Do not leave it in contact with dampness for prolonged periods. Follow these steps for cleaning: 1st step: Remove drill from contra-angle (if applicable). 2nd step: Immerse drill completely in enzymatic detergent (10% solution). 3rd step: Wash it in an ultrasonic washer for approximately 10 minutes. 4th step: Rinse it with plenty of distilled water, until all solution residues have been completely removed. The use of nylon brushes is permitted. 5th step: Dry it with a paper towel or compressed air. 6th step: Perform a visual inspection, observing whether there are any faults in the cleaning process. If you still find any residues, reimmerse the drill in detergent - 2nd step - and clean it using a nylon brush, if necessary. Repeat rinsing and drying sequence. 7th step: If applicable, place the drill in an appropriate, properly cleaned and dried case. 8th step: Select packaging according to sterilization process. Preferably use self-sealing package of medical grade paper with laminated film. IMPORTANT: Never store drills without drying them thoroughly. Drying is very important to prevent oxidation. CAUTION: Do not use incrustation removal products, as they tend to darken drills and cause oxidation. Avoid using enzymatic liquid in concentrations above 10%, and drying drills that still contain residues from cleaning solutions, as these procedures favor oxidation. The use of steel brushes is also not recommended.

The drying of the parts is of utmost importance before storage and sterilization, because the accumulation of moisture on the products is harmful and may cause oxidation.

**PRESENTATION AND STERILIZATION**

These products are reusable and supplied non-sterile, packed individually in a package that provides double protection. These products must be correctly sanitized and sterilized before use. Sterilize the products the day before or on the day of the procedure. ATTENTION: It is not recommended to autoclave these products in their original packaging.

We preferably recommend following the sterilization method by steam autoclaving and, yet, the parameters and procedures established in standard BS EN ISO 17665-1:2006 (Sterilization of health care products. Moist heat Requirements for the development, validation and routine control of a sterilization process for medical devices). Alternatively, can proceed according to the instruction manual of the manufacturer of the autoclave. Suitable parameters: 121 °C, 30 min, at 1 atm. Sterilization validity: 7 to 15 days, provided that the product is conditioned in clean, dry and protected from the sunlight environment.

## PRECAUTIONS

- Inappropriate planning can jeopardize the performance of the implant unit leading to failures of the system, as loss or fracturing of the implant.
- The Guided Surgery technique allows for the installation of Titamax Morse Taper, Alvim Morse Taper, Drive Morse Taper, and Facility Implants.
- The Guided Surgery technique does not allow for installation of implants with a diameter equal to or greater than 5.0 mm or with a length exceeding 17 mm.
- To perform the procedure, the patient must have sufficient interocclusal space in the area for installation of the implant, for manipulation of the instruments of the Guided Surgery technique.
- Be aware in cases patients that show signs of allergy or hypersensitivity to the chemical components of the material: surgical stainless steel
- Do not use the product if the packaging is broken.
- Do not use the product with the validity expired.
- Before each procedure, make sure the pieces are properly seated.
- Be sure you have all the instruments needed to perform the procedure according to the surgical planning.
- Before each procedure check the conditions of the Neodent surgical instruments, always respecting their useful life. Replace the instruments if there is damage, markings deleted,

sharpening jeopardized, deformation and wear.

- During the surgical procedure, make sure you don't exert lateral movements with the surgical instruments, avoiding the loss of the mounting of the Sleeve.
- Ensure that the parts are not swallowed or aspirated by the patient.
- Always use the sequence of Neodent products for the Neodent Guided Surgery technique. The use of instruments or Sleeves of other manufacturers does not ensure the perfect function of the Neodent Implant System and exempts any product warranty.
- Drills should not be resharpened.
- It is the professional's responsibility to use the Neodent products according to the instructions for use.

## ADVERSE EFFECTS

No adverse effects will occur, provided that the indications are complied with.

## POST-OPERATIVE PRECAUTIONS AND MAINTENANCE

Instruct the patient as to the need of a professional medical monitoring after the surgery and to obey the guidelines regarding the precautions, hygiene and prescription of drugs. These guidelines are the responsibility of the professional in charge.

## STORAGE CONDITIONS

This product should be stored in a clean and dry location, in a maximum temperature of 40 °C and

protected from direct sunlight.

## DISPOSAL OF MATERIAL

After use, all the products and consumables used in the surgery for the installation of dental implants may endanger the health of anyone who handles them. Before discarding them into the environment, it is recommended to take a look at the current legislation and adhere to it.

## DATE OF EXPIRATION

Written on the label.

## LIFE CYCLE

The Guided Surgery Drill Guide must be discarded when it loses its functionality. Guided Surgery Drills and Guided Surgery Pilot Drills have a useful life of up to 20 drillings for bone qualities I and II and up to 25 drillings for bone qualities III and IV, provided that its conditions for use are respected, including proper irrigation and rotation. Regardless of the lifespan mentioned, the professional should always check the sharpness condition of the drill. Drills with improper conditions for cutting should be discarded. Note: Classification of the bone quality according to Lekholm and Zarb (1985).

Federal (USA) law restricts this device to sale by or on the order of a licensed dentist or physician.

Not all products are available in all countries. Please, contact the authorized distributor.

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#### Titamax CM EX

	Guided Surgery 2.0	Guided Surgery 2.8	Guided Surgery 3.0
3.5 mm	✓		
3.75 mm	✓	*	
4.0 mm	✓	*	*

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\* Opcional

#### Titamax CM

	Broca Helicoidal Plus 2.0	Broca Helicoidal Plus 2.8	Broca Helicoidal Plus 3.0	Broca Piloto Guided Surgery 3.5	Broca Helicoidal Plus 3.3	Broca Piloto Guided Surgery 3.75	Broca Piloto Guided Surgery Plus 4.0
3.5 mm	✓	✓		✓			
3.75 mm	✓		✓			✓	
4.0 mm	✓		✓		✓		✓

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#### Drive CM

	Broca Guided Surgery Alvim 2.0	Broca Guided Surgery Alvim 3.5	Broca Piloto Guided Surgery 3.5	Broca Guided Surgery Alvim 4.3	Broca Piloto Guided Surgery 34.3
3.5 mm	✓	✓	*		
4.3 mm	✓	✓		✓	*

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\* Opcional

#### Alvim CM

	Broca Guided Surgery Alvim 2.0	Broca Guided Surgery Alvim 3.5	Broca Piloto Guided Surgery 3.5	Broca Guided Surgery Alvim 4.3	Broca Piloto Guided Surgery 4.3	Macho de Rosca Guided Surgery Alvim 3.5	Macho de Rosca Guided Surgery Alvim 4.3
3.5 mm	✓	✓	*			**	
4.3 mm	✓	✓		✓	*		**

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\*\* Recomendada en tipos óseos I y II

#### Facility

	Broca Guided Surgery Facility 2.0	Broca Guided Surgery Facility 2.0/2.4	Broca Piloto Guided Surgery Facility	Macho Guided Surgery Facility
10 mm	✓	✓	*	**
12 mm	✓	✓	*	**
14 mm	✓	✓	*	**

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\* Opzionale

Titamax CM

	Fresa a Spirale 2.0	Fresa a Spirale 2.8	Fresa a Spirale 3.0	Fresa Pilota Guided Surgery 3.5	Fresa a Spirale 3.3	Fresa Pilota Guided Surgery 3.75	Fresa Pilota Guided Surgery 4.0
3.5 mm	✓	✓		✓			
3.75 mm	✓		✓			✓	
4.0 mm	✓		✓		✓		✓

✓ Sequenza consigliata

Drive CM

	Fresa Guided Surgery Alvim 2.0	Fresa Guided Surgery Alvim 3.5	Fresa Pilota Guided Surgery Alvim 3.5	Fresa Guided Surgery Alvim 4.3	Fresa Pilota Guided Surgery Alvim 4.3
3.5 mm	✓	✓	*		
4.3 mm	✓	✓		✓	*

✓ Sequenza consigliata

\* Opzionale

Alvim CM

	Fresa Guided Surgery Alvim 2.0	Fresa Guided Surgery Alvim 3.5	Fresa Pilota Guided Surgery Alvim 3.5	Fresa Guided Surgery Alvim 4.3	Fresa Guided Surgery Alvim 4.3	Maschiatore Alvim Guided Surgery 3.5	Maschiatore Alvim Guided Surgery 4.3
3.5 mm	✓	✓	*			**	
4.3 mm	✓	✓		✓	*		**

✓ Sequenza consigliata

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Facility

	Fresa Guided Surgery Facility 2.0	Fresa Guided Surgery Facility 2.0/2.4	Fresa Pilota Guided Surgery Facility	Maschiatore Facility Guided Surgery
10 mm	✓	✓	*	**
12 mm	✓	✓	*	**
14 mm	✓	✓	*	**

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\* En option

Titamax CM

	Foret Guided Surgery 2.0	Foret Guided Surgery 2.8	Foret Guided Surgery 3.0	Foret Pilote Guided Surgery 3.5	Foret Guided Surgery 3.3	Foret Pilote Guided Surgery 3.75	Foret Pilote Guided Surgery 4.0
3.5 mm	✓	✓		✓			
3.75 mm	✓		✓			✓	
4.0 mm	✓		✓		✓		✓

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Drive CM

	Foret Guided Surgery Alvim 2.0	Foret Guided Surgery Alvim 3.5	Foret Pilote Guided Surgery 3.5	Foret Guided Surgery Alvim 4.3	Foret Pilote Guided Surgery 4.3
3.5 mm	✓	✓	*		
4.3 mm	✓	✓		✓	*

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\* En option

Alvim CM

	Foret Guided Surgery Alvim 2.0	Foret Guided Surgery Alvim 3.5	Foret Pilote Guided Surgery 3.5	Foret Guided Surgery Alvim 4.3	Foret Pilote Guided Surgery 4.3	Taraud pour os Guided Surgery Alvim 3.5	Taraud pour os Guided Surgery Alvim 4.3
3.5 mm	✓	✓	*			**	
4.3 mm	✓	✓		✓	*		**

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Facility

	Foret Guided Surgery Facility 2.0	Foret Guided Surgery Facility 2.0/2.4	Foret Pilote Guided Surgery Facility	Taraud pour os Facility Guided Surgery
10 mm	✓	✓	*	**
12 mm	✓	✓	*	**
14 mm	✓	✓	*	**

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