ST.330.326

Português

Transfer de Moldeira Aberta e Fechada

Italiano

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This device is intended for a specialized procedure, which should be performed by professionals qualified in Dental Implants. For optimum results, use the product knowing the appropriate techniques. Always apply them under appropriate conditions.

DESCRIPTION

The Implant Impression Coping comes with a screw and is made of titanium alloy, in accordance with ASTM standard F136. It has a cylindrical shape with concave facets. It is available in the regular or long models, with CM interface.

The Implant Impression Coping for the closed tray technique comes with a pair of Positioner Caps made of POM (Polyacetal).

INDICATIONS FOR USE

This product is indicated for transferring the position and orientation of the implant installed in the mouth.

APPLICATIONS

This produce allows you to transfer, via impression, the three-dimensional position of the implant installed in the mouth. The solution is presented for impression techniques with a closed or open tray.

The Implant Impression Coping is indicated according to the table below:

Prosthetic interface	Form of presentation
CM	Rotational
CM Exact	Anti-rotational

CONTRAINDICATIONS

The CM Implant Exact Impression Coping models are not indicated for multiple restorations. This product is contraindicated for patients exhibiting signs of allergy or hypersensitivity to the chemical ingredients of the material: Titanium alloy Ti6Al4V-ELI and POM.

HANDLING

The Implant Impression Coping must be selected according to prior prosthetic planning, respecting the clinical situation of the patient.

Open Tray: fit the body of the Impression Coping together with the implant selected, and turn only the screw, manually or with the aid of the 1.2 Connection / 1.2 Wrench until the whole assembly is stable. Do not exceed the maximum torque of 10 N.cm). For multiple restorations, the joint of the impression copings must be done with acrylic resin. Make an opening in the

impression tray to allow for coupling of the 1.2 Connection / 1.2 Wrench after the impression is taken. Handle the impression material according to the recommendations of the manufacturer of the material. Load the impression tray with the material and bring it to the patient's mouth. After the material sets, use the 1.2 Connection / 1.2 Wrench to loosen the screw. Remove the impression tray and make sure the Impression Coping has come with the mold. Carefully fit and tighten the analog of the implant onto the Impression Coping. Make sure you do not move the Impression Coping during the fitting of the analog. Create the plaster model according to the proper techniques. Remove the Transfer from the plaster model.

Closed Tray: fit the body of the Impression Coping together with the implant selected, and turn only the screw with the aid of the 1.2 Connection / 1.2 Wrench until the whole assembly is stable. Do not exceed the maximum torque of 10 N.cm). Fit the Positioner Cap over the Impression Coping. Make sure that it is well seated. Handle the impression material according to the recommendations of the manufacturer of the material. Load the impression tray with the material and bring it to the patient's mouth. After the material sets, remove the impression tray. The Positioner Cap should come with the mold. With the aid of the 1.2 Connection / 1.2 Wrench, loosen the screw and remove the Implant Impression Coping. Then screw the Impression Coping onto the analog of the implant. Carefully fit the coping-analog assembly in the mold. Ensure a correct fit between the Impression Coping and the Positioner Cap. Create the plaster model according to the proper techniques. Remove the Transfer from the plaster model.

PRESENTATION AND STERILIZATION

This product is indicated for single use and provided in non-sterile form, and must therefore be sterilized before use. It is packed in a blister-pack containing one impression coping with screw and two Positioner Caps.

For sterilization, all the items must be removed from their original packaging.

During sterilization, the Positioner Caps should not be held in the same package as the impression coping and screw.

It is recommended, preferably, to follow the method of sterilization by steam autoclave (moist heat), gravitational or vacuum autoclave, with exposure time of 3 minutes at 132° C (270 ° F). The product must be accommodated on an appropriate tray. Use it immediately after sterilization; do not store.

PRECAUTIONS

- The Open Tray technique is not recommended for patients with narrow opening of the mouth.
- The CM prosthetic interface is not compatible with the WS prosthetic interface.
- Make sure that the stability of the implant is sufficient to support the procedures of insertion and removal of the Implant Impression Coping.
- Surgical planning and/or inadequate prosthesis can compromise the performance of the implant/prosthesis assembly, resulting in system failure, such as loss or fracture of the implant, loosening or fracture of components and/or prosthetic screws.
- This product is of single use and cannot be reprocessed.
- Reuse of this product may cause: adverse biological effects arising from product residue, microorganisms, and/or substances originating from previous uses and/or reprocessing; changes to the original macro• and micro-structural physical, mechanical, and chemical characteristics of the product, which may affect their intended functionality;
- Reuse of this product does not guarantee its safety and effectiveness, and invalidates any warranty for related products.
- Do not use the product if the packaging is broken.
- Do not use the product after expiry date.
- For immediate load application, check the torque indication of the implant installed.
- See the torque to be given on the prosthetic component to be used. The excess or lack of torque can lead to undesirable result.
- Before each procedure, make sure the pieces are properly seated.
- Ensure that the parts are not swallowed or aspirated by the patient.
- Before each procedure check the conditions of the Neodent surgical instruments, always respecting their useful life. Replace the instruments if there is damage, markings deleted, sharpening jeopardized, deformation and wear.
- Always use the Neodent product sequence. The use of prosthetic components and/or instruments of other manufacturers does not ensure the perfect function of the Neodent Implant System and exempts any product warranty.
- It is the professional's responsibility to use the Neodent products according to the instructions for use.

ADVERSE EFFECTS

No adverse effects will occur, provided that the indications are complied with.

POST-OPERATIVE PRECAUTIONS AND MAINTENANCE

Instruct the patient as to the need of a professional medical monitoring after the procedure and to obey the guidelines regarding the precautions, hygiene and prescription of drugs. These guidelines are the responsibility of the professional in charge.

STORAGE CONDITIONS

This product should be stored, in its original packaging, in a clean and dry location, in a maximum temperature of 40°C and protected from direct sunlight.

DISPOSAL OF MATERIAL

After use, all the products and consumables used in the surgery for the installation of dental implants may endanger the health of anyone who handles them. Before discarding them into the environment, it is recommended to take a look at the current legislation and adhere to it.

DATE OF EXPIRATION

Written on the label.

Federal (USA) law restricts this device to sale by or on the order of a licensed dentist or physician.

Not all products are available in all countries. Please, contact the authorized distributor.

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