English	MRI Safety Information (for Europe)	2 – 3
Deutsch	MRT-Sicherheitsinformation (für Europa)	4 – 5
Français	Sécurité dans l'environnement de résonance magnétique (pour l'Europe)	6 – 7
Italiano	Informazioni di sicurezza relative alla risonanza magnetica per immagini (MRI) (per l'Europa)	8 – 9
Español	Información de seguridad para resonancias magnéticas (para Europa)	10 – 11
Português	Informações sobre segurança em ambiente de RM (para a Europa)	12 – 13

Magnetic Resonance Imaging (MRI) has become a widely used and accepted tool for medical imaging. It uses magnetic fields and radio waves to generate images of the body. As dental implants and prostheses are made of or include metal parts, concerns regarding the safety of such appliances in the high magnetic fields used are frequently raised by patients and dental practitioners. Main concerns include magnetically induced displacement forces, torque and radio frequency (RF) induced heating capable of compromising patient safety and image artifacts, potentially interfering with image quality.

Straumann titanium implants are made of titanium (Ti) quality grade 4, or Roxolid®, a titanium-zirconium (TiZr) alloy containing around 85% titanium. These materials are considered paramagnetic and therefore only weakly interact with magnetic fields.¹

Straumann® secondary components, such as abutments, copings, closure screws and healing caps are made from titanium grade 2 and 4, TAN (titanium-6-aluminium-7-niobium), TAV (titanium-6-aluminium-4-vanadium), gold-based alloys like Ceramicor®, Elitor®, Esteticor® (registered trademarks of Cendres and Métaux SA, Biel; Switzerland), PEEK, PMMA or ceramics (ZrO₂). These materials are either paramagnetic or non-magnetic. Coron®, a cobalt-chromium alloy frequently used for the production of prosthetics is generally considered as non-magnetic.

Straumann has conducted a series of studies to assess the behavior of the devices of the Straumann® Dental Implant System when subjected to the magnetic fields of MRI equipment. Following is a summary of the results of these studies. This is provided as a guide for the radiologist in assessing whether patients fitted with Straumann implants and restorations should be subjected to MRI scans.

MRI Safety evaluations of various materials on test specimens (magnetic displacement and artifacts)

This study assessed displacement force effects in the main static magnetic field of a 3 Tesla (T) MRI scanner per ASTM F2052-14 as well as calculating the expected magnetic force in a stronger magnetic field gradient of 300 mT/cm per ASTM F2213-06 (2011)².

Assessment of image artifacts at 3 T according to ASTM F2119-07 (2013) was also conducted. The test samples consisted of $\varnothing 5 \times 10$ mm cylinders of the following materials:

Group	Materials
Non-conductive	PEEK polymer, PMMA (Vita CAD Temp, Polycon ae)
Gold alloys	Ceramicor®, Ney H, Alloy 6019
Ti and Ti Alloys	Titanium (grade 2 and 4), Ti-13Zr (Roxolid®), Ti-6Al-7Nb, Ti-6Al-4V
CoCr	CoCr alloy

The tested samples all showed magnetic displacement forces less than gravitational forces. Deflection angles were $\leq 30^\circ$. The magnetic displacement forces of less than 10 mN (weight of worst case sample: 1.72 g) do not pose any safety risk. No torque effects could be seen for the tested samples in a magnetic field up to 3T. The magnitude of image artifacts is based upon the implant material.

We identified four different material groups exhibiting image artifacts of similar magnitude:

Group	Materials	Gradient Echo Artifact Magnitude	Spin Echo Artifact Magnitude
Non-conductive	PEEK polymer, PMMA	1 mm	1 mm
Gold alloys	Ceramicor, Ney H, Alloy 6019	7 to 10 mm	5 to 10 mm
Ti and Ti Alloys	Titanium, Ti-13Zr (Roxolid®), Ti-6Al-7Nb, Ti-6Al-4V	14 to 16 mm	12 to 15 mm
CoCr	CoCr alloy	22 to 24 mm	20 to 24 mm

MRI Safety Evaluations of dental implant constructions (magnetic displacement, RF heating and artifacts)

This study assessed displacement force effects in the main static magnetic field of a 3 T MRI scanner per ASTM F2052-14 as well as calculating the expected magnetic force in a stronger magnetic field gradient of 300 mT/cm per ASTM F2213-06 (2011)³.

Assessment of image artifacts at 3 T according to ASTM F2119-07 (2013) was also conducted; and heating effects due to the radio frequency field (RF-field) during MRI scans at 1.5 T and 3 T according to ASTM F2182-11a were assessed.

The test samples consisted of the following:

- A sixteen unit titanium grade 4 bridge mounted to sixteen implants constructed of Ti-13Zr alloy
- A sixteen unit CoCr alloy bridge mounted to sixteen implants constructed of Ti-13Zr alloy
- A full arch titanium grade 4 bar assembly mounted to four implants constructed of Ti-13Zr alloy
- A full arch CoCr alloy bar assembly mounted to four implants constructed of Ti-13Zr alloy

As seen in the first study, the tested samples all showed deflection angles of $\leq 30^\circ$. Magnetic displacement forces were less than 100 mN and do not pose any safety risk. No torque effects could be seen for the tested samples in a magnetic field up to 3 T.

The titanium samples showed artifacts up to 26 mm from the surface, whereas the CoCr samples showed artifacts up to 42 mm. Spin echo sequences showed reduced artifacts of up to 13 mm (Ti) and up to 20 mm (CoCr).

RF heating tests in worst case orientation as described by the ASTM standard and with a whole body specific absorption rate (SAR) of 2 W/kg (normal operating mode) showed temperature increases of up to 5.2 °C (1.5 T) and up to 5.9 °C (3 T) for the tested implants. No cooling effects from perfusion and blood flow have been taken into account.

Summary of the non-clinical testing:

Magnetic displacement:

- The Straumann construct exhibited displacement forces of less than 100 mN (in a static magnetic field of 3 T and a whole body averaged SAR of 2 W/kg in normal operating mode).

RF heating:

- The metallic components of the Straumann® Dental Implant System are expected to produce a maximum temperature rise of less than 5.9 °C at 3 T and of less than 5.2 °C at 1.5 T (after 15 min continuous scanning).

Image artifacts:

In a 3 T MRI system, the approximate extension of image artifacts for full arch assemblies was:

- 26 mm (Ti) or 42 mm (CoCr) when imaged using gradient echo
- 13 mm (Ti) or 20 mm (CoCr) when imaged using spin echo respectively

Conclusion

Based on non-clinical testing it can be concluded that the components of the Straumann® Dental Implant System are unlikely to interfere with patient safety under the MRI conditions defined above. Magnetic displacement of components of the Straumann® Dental Implant System has been shown to be less than the force exerted on the device by gravity, and RF heating leads to a maximum temperature increase of less than 6 °C, not taking into account the cooling effect of surrounding tissue and blood flow. However, image artifacts are to be expected and have to be considered at image analysis.

Disclaimer

Note that these studies have not been found sufficient to assign a rating of MR Safe, MR Conditional or MR Unsafe to the components of the Straumann® Dental Implant System; scanning a patient fitted with these devices may result in patient injury.

Due to the large variety of MRI scanners available on the market, Straumann cannot make any predictions regarding the safety or behavior of our implants and components in any specific MRI system.

Finally, Straumann cannot take any responsibility for the composition and behavior of any third party product (including crown, bridge, bar, denture, etc.), which is not distributed by Straumann and may contain materials which may not be compatible with MRI imaging.

¹ P. A. Tipler, "Physics for Scientists and Engineers", Sixth edition, Publisher: Palgrave Macmillan, 2008, ISBN-13: 9781429202657

² Test Report: MRI Safety Evaluations of Various materials. Report number: GT_TR_straumann_2014_01 Rev A, 22 January 2015

³ Test Report: MRI Safety Evaluations of Dental Implants. Report number: GT_TR_straumann_2014_02 Rev A, 23 January 2015

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