ST.330.274

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This device is intended for a specialized procedure, which should be performed by professionals qualified in Dental Implants. For optimum results, use the product knowing the appropriate techniques. Always apply them under appropriate conditions, also in an operating room atmosphere.

DESCRIPTION

Instrument manufactured in stainless steel. In its active end it has a tip supported by Facility Prosthetic Components. At the other end, it has a spring system activated axle.

APPLICATIONS


It is suitable for installation of the Prosthetic Components Facility line by means of impact.

CONTRAINDICATIONS

This product does not show contraindications, provided that it be properly used for the shown purposes.

HANDLING

For disassembly prior to cleaning, unscrew the Placement Aid's body, releasing all its internal parts. Check the figure:

For sterilization, reassemble Placement Aid's parts.

1. Attach the spring (2) to the axle (3) and insert this set to the disassembled part of the greater length body (1). 2. Attach the tip (4) to the other part of the body (shorter length)(5). 3. Screw body parts. 4. Ensure Facility Abutment Placement Aid functioning by pulling the axle (3). With the Facility Prosthetic Component already on the Implant, place the tip (4) of the Placement Aid on it. The Placement Aid must be aligned along the axis of the Prosthetic Component. Keep Placement Aid slightly pressed, ensuring contact between the tip and the Prosthetic Component during impact. Pull the axle until the end of the spring stroke and then release it. We recommend 3 impacts for an effective component installation.

SANITATION

This product must be properly cleaned after each use. For cleaning, do the following: 1st step: Disassemble the Facility Abutment Placement Aid by unscrewing its body and releasing all parts. Keep the O'ring positioned in Placement

Aid body's thread. 2nd step: Fully immerse the part in enzymatic detergent (diluted according to the manufacturer) 3rd step: Wash it in an ultrasonic cleaner for approximately 10 to 15 minutes. 4th step: Rinse with distilled water in abundance, until completely remove the wastes from the solution. Nylon brushes use is allowed. 5th step: Dry it with a clean, dry cloth or with compressed air. 6th step: Perform a visual inspection, observing if there are failures in the cleaning process. If there are still wastes, the part should be once again immersed in detergent - 2nd step - and, if necessary, the cleaning should be performed with the aid of a nylon brush. Repeat the rinsing and drying sequence. 7th step : Select the package according to the sterilization process. Use preferably a self-sealing package if surgical grade paper with laminated film. ATTENTION: Do not use anticrust products, avoid using enzymatic liquid in concentrations above 10% and dry parts that still have cleaning solution wastes, because these procedures cause oxidation. The use of steel brushes is also inadvisable.

PRESENTATION AND STERILIZATION

This product is reusable and supplied unitarily, non-sterile, conditioned in envelope (film type) packing and, therefore, should be sterilized before use. It must be sterilized after being mounted. Sterilize the product on the day before or on the day of the procedure. We preferably recommend following the sterilization method by steam autoclaving and, yet, the parameters and procedures established in standard BS EN ISO 17665-1:2006 (Sterilization of health care products. Moist heat. Requirements for the development, validation and routine control of a sterilization process for medical devices). Alternatively, can proceed according to the instruction manual of the manufacturer of the autoclave. Suitable parameters: 121 °C, 30 min, at 1 atm. Sterilization validity: 7 to 15 days, provided that the product is conditioned in clean, dry and protected from the sunlight environment.

PRECAUTIONS

- Do not use the product if the package has been violated.
- This product should be sterile.
- This product should be disassembled for cleaned and assembled for sterilization.
- Ensure Placement Aid functioning by pulling the axle.
- Autoclaving the product in its original packaging is not recommended.
- This product

should be handled firmly, avoiding contact with adjacent soft tissues.

- Do not remove the O'ring from Placement Aid's body to sanitize or sterilize.
- Better results are achieved with the utilization of the Neodent products sequence. The utilization of instrumentations and/or Prosthetic Components from other manufacturers does not assure the perfect function of the Neodent Implant System and voids any product guarantee.
- Improper planning can jeopardize the performance of the Implant/Prosthesis set, causing failures to the system, such as loss or break of the Implant, loosening or fracture of the Prosthetic Screws.

ADVERSE EFFECTS

Adverse effects will only occur if the choice and the use of instrument are improper.

POST-OPERATIVE PRECAUTIONS AND MAINTENANCE

Instruct the patient as to the need of a professional medical monitoring after the surgery and to obey the guidelines regarding the precautions, hygiene and prescription of drugs. These guidelines are the responsibility of the professional in charge.

STORAGE CONDITIONS

This product should be stored in a clean and dry location, in a maximum temperature of 40°C and protected from direct sunlight.

DISPOSAL OF MATERIAL

Every product and consumable used during the surgery for the installation of dental implants may endanger the health of those who handle them after use. Before discarding them into the environment, it is recommended to take a look at the current legislation and adhere to it.

DATE OF EXPIRATION

Written on the label.

LIFE CYCLE

This product should be disposed of when it loses its functionality.

Federal (USA) law restricts this device to sale by or on the order of a licensed dentist or physician.

Not all products are available in all countries. Please, contact the authorized distributor.

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